

DIAMOND GRIP PLUS™

POWDER-FREE LATEX EXAMINATION GLOVES

● STANDARD EXAMINATION GLOVE



Features & Benefits:

- Polymer-coated for easy, powder-free donning
- Longer and thicker than typical exam gloves for added protection
- Fully textured for a secure grip in wet or dry conditions
- Premium, low-modulus natural rubber latex for a true fit and reduced hand fatigue
- Reduced protein
- EN 374-2, EN 374-3 certified



Performance and comfort is a premium.

Diamond Grip Plus™ is made from soft, premium natural rubber latex for a comfortable fit. Fully textured for enhanced sensitivity and control. Manufactured with integrated polymer technologies for powder-free donning and product consistency. **Diamond Grip Plus™** meets FDA and EN/ISO standards for exam glove quality and performance.

EN 374



MICROFLEX®
THE MOST TRUSTED NAME IN GLOVES®

www.microflex.com/int

REGISTERED
ISO 13485

DIAMOND GRIP PLUS™

POWDER-FREE LATEX EXAMINATION GLOVES

STANDARD EXAMINATION GLOVE



GLOVE PROFILE

Diamond Grip Plus*

Length	248mm / 9.8 in.
Cuff Thickness	3.9 mils
Palm Thickness	5.1 mils
Finger Thickness	6.3 mils
Force at Break	min. 9
Elasticity/Elongation (%)	min. 700
Powder Mass	max. 2 mg / glove
Protein Level	max. 50 µg/g

GLOVE ATTRIBUTES

Material:	Premium quality natural rubber latex
Type:	Non-sterile
Shape:	Ambidextrous
Color:	Natural
Sizes:	Extra-Small, Small, Medium, Large, Extra-Large
Cuff:	Beaded
Glove Interior:	Polymer-coated, powder-free
Glove Exterior:	Fully textured
Usage:	For single use only
Packaging:	100 gloves per dispenser box, by weight; 10 boxes per case

* ISO1193-1; EN455 I, II, III; EN 374-2, EN 374-3; CE Class I

Powder-free latex examination gloves

See for yourself why **Microflex®** gloves are the first choice of exam glove users nationwide. For samples or more information on any of the fine latex and synthetic products from **Microflex®**, contact your **Microflex®** distributor or for the name of a distributor in your area, call:

775-746-6600 ext. 1726

Or visit our website at:

www.microflex.com/int

Fax: 775-746-6586.

PRODUCT CODE	RANGE OF SIZE	DESCRIPTION	CASE QUANTITY
DGP-INT-XS	EXTRA-SMALL	Latex, Powder-free, Examination	1000
DGP-INT-S	SMALL	Latex, Powder-free, Examination	1000
DGP-INT-M	MEDIUM	Latex, Powder-free, Examination	1000
DGP-INT-L	LARGE	Latex, Powder-free, Examination	1000
DGP-INT-XL	EXTRA-LARGE	Latex, Powder-free, Examination	1000

Caution Statement: This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.



MICROFLEX®
THE MOST TRUSTED NAME IN GLOVES®

www.microflex.com/int

©2007 Microflex Corporation
P.O.Box 32000
Reno, NV 89533-2000
All Rights Reserved.
07.00003 AL 061907