


# MICRO ONE®

## LIGHTLY POWDERED LATEX EXAMINATION GLOVES

### • STANDARD EXAMINATION GLOVE

#### Features & Benefits:

- Thicker and longer than typical examination gloves for added protection where you need it most
- Smooth surface for excellent tactile sensitivity
- Lightly powdered for easy donning
- Meets EN/ISO requirements for quality and performance
- EN 374-2, EN 374-3 certified



**Lightly powdered.  
A proven performer.**

**Micro One®** is a highly versatile examination glove that performs well in demanding conditions. Made from soft latex, **Micro One®** has a natural fit for reduced hand fatigue. Lightly powdered to go on easily with minimal excess powder. Smooth surface allows for maximum sensitivity. Thicker and longer than typical latex exam gloves for added protection where you need it most.

EN 374



**MICROFLEX®**  
THE MOST TRUSTED NAME IN GLOVES®

[www.microflex.com/int](http://www.microflex.com/int)

REGISTERED  
**ISO 13485**

# MICRO ONE®

## LIGHTLY POWDERED LATEX EXAMINATION GLOVES

### ● STANDARD EXAMINATION GLOVE

#### GLOVE PROFILE

#### Micro One\*

Length	245mm / 9.6 in.
Cuff Thickness	3.9 mils
Palm Thickness	5.5 mils
Finger Thickness	6.3 mils
Force at Break	min. 9
Elasticity/Elongation (%)	min. 750
Powder Mass	max. 10 mg / dm <sup>2</sup>
Protein Level	max. 150 µg/g

#### GLOVE ATTRIBUTES

Material:	Premium quality natural rubber latex
Type:	Non-sterile
Shape:	Ambidextrous
Color:	White
Sizes:	Extra Small, Small, Medium, Large, Extra Large
Cuff:	Beaded
Glove Interior:	Lightly powdered
Glove Exterior:	Smooth, non-textured
Usage:	For single use only
Packaging:	100 gloves per dispenser box, by weight; 10 boxes per case

\* ISO11193-1; EN455 I, II, III; EN 374-2, EN 374-3; CE Class I



### Powder-free latex examination gloves

See for yourself why **Microflex®** gloves are the first choice of exam glove users nationwide. For samples or more information on any of the fine latex and synthetic products from **Microflex®**, contact your **Microflex®** distributor or for the name of a distributor in your area, call:

**775-746-6600** ext. 1726

Or visit our website at:

[www.microflex.com/int](http://www.microflex.com/int)

Fax: 775-746-6586

PRODUCT CODE	RANGE OF SIZE	DESCRIPTION	CASE QUANTITY
MO-INT-XS	EXTRA-SMALL	Latex, Lightly Powdered, Examination	1000
MO-INT-S	SMALL	Latex, Lightly Powdered, Examination	1000
MO-INT-M	MEDIUM	Latex, Lightly Powdered, Examination	1000
MO-INT-L	LARGE	Latex, Lightly Powdered, Examination	1000
MO-INT-XL	EXTRA-LARGE	Latex, Lightly Powdered, Examination	1000

**CAUTION STATEMENT:** This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.



**MICROFLEX®**  
THE MOST TRUSTED NAME IN GLOVES®

[www.microflex.com/int](http://www.microflex.com/int)

©2007 Microflex Corporation  
P.O.Box 32000  
Reno, NV 89533-2000  
All Rights Reserved.  
07.00003 AL 071207