

# SUPRENO PLUS™


POWDER-FREE NITRILE EXAMINATION GLOVES

• FOR HIGH RISK ENVIRONMENTS



## Features & Benefits:

- Soft and durable nitrile formulation for improved elasticity and greater comfort
- Interior polymer coating makes glove easy to don
- 100% nitrile to reduce the potential for adverse reactions associated with latex protein
- Textured fingers for reliable grip in wet and dry conditions
- Extended cuff and gradient dipping provides added protection where you need it most
- Meets all current FDA and EN/ISO standards assuring consistent quality in every box, every time
- EN 374-2, EN 374-3 certified



**Exceptional comfort in a tough nitrile glove.**

**Supreno Plus™** is a soft and durable nitrile examination glove for high-risk applications. This strong nitrile formulation makes for a comfortable, true-fitting glove that won't fatigue hands over long periods of wear. An extended cuff protects the wrist and forearm, and textured fingers offer better wet and dry gripping. Interior polymer coating allows **Supreno Plus™** to don easily. Online chlorination processing means consistent protection in every box, every time.

Test data available upon request



BS EN 374-2  
Micro-organism  
Hazards



BS EN 374-3  
Chemical  
Hazards

**MICROFLEX®**  
THE MOST TRUSTED NAME IN GLOVES®  
[www.microflex.com/int](http://www.microflex.com/int)

REGISTERED  
**ISO 13485**

# SUPRENO PLUS™

## POWDER-FREE NITRILE EXAMINATION GLOVES

### FOR HIGH RISK ENVIRONMENTS



#### GLOVE PROFILE

##### Supreno Plus\*

Length	295mm / 11.6 in.
Cuff Thickness	4.7 mils
Palm Thickness	5.9 mils
Finger Thickness	7.9 mils
Force at Break	min. 9
Elasticity/Elongation (%)	min. 500
Powder Mass	max. 2 mg / glove

#### GLOVE ATTRIBUTES

<b>Material:</b>	Nitrile; not made from natural rubber latex
<b>Type:</b>	Non-sterile
<b>Shape:</b>	Ambidextrous
<b>Color:</b>	Blue
<b>Sizes:</b>	Small, Medium, Large, Extra-Large
<b>Cuff:</b>	Extended, beaded
<b>Glove Interior:</b>	Polymer-coated, powder-free
<b>Glove Exterior:</b>	Textured fingers
<b>Usage:</b>	For single use only
<b>Packaging:</b>	50 gloves per dispenser box, by weight; 10 boxes per case

\* ISO1193-1; EN455 I, II, III; EN420; EN374-2; EN374-3, CE Class I

### Powder-free synthetic examination gloves

See for yourself why **Microflex®** gloves are the first choice of exam glove users nationwide. For samples or more information on any of the fine latex and synthetic products from **Microflex®**, contact your **Microflex®** distributor or for the name of a distributor in your area, call:

**775-746-6600** ext. 1726

Or visit our website at:

[www.microflex.com/int](http://www.microflex.com/int)

Fax: 775-746-6586

PRODUCT CODE	RANGE OF SIZE	DESCRIPTION	CASE QUANTITY
SP-INT-S	SMALL	Nitrile, Non-sterile, Examination	500
SP-INT-M	MEDIUM	Nitrile, Non-sterile, Examination	500
SP-INT-L	LARGE	Nitrile, Non-sterile, Examination	500
SP-INT-XL	EXTRA-LARGE	Nitrile, Non-sterile, Examination	500

**Caution Statement:** Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.



**MICROFLEX®**  
THE MOST TRUSTED NAME IN GLOVES®

[www.microflex.com/int](http://www.microflex.com/int)

©2007 Microflex Corporation  
P.O. Box 32000  
Reno, NV 89533-2000  
All Rights Reserved.  
07.00003 AL 061907

EUROPEAN PATENT NO. 0456333.