

SYNETRON® PLUS

**POWDER-FREE LATEX
EXTENDED CUFF EXAMINATION GLOVES**

Recommended for use in the following industries:


EMS / Laboratory / Industrial Safety

Synetron® Plus latex examination gloves are a strong barrier against many chemicals and biohazards, and its comfortable during extended periods of wear. Synetron® Plus is textured for a secure grip and is polymer-coated for product consistency and easy donning.

FEATURES AND BENEFITS:

- Interior/exterior polymer-coating system provides superior anti-sticking properties that assist with highly effective donning in wet, damp or dry conditions
- Extended cuff for increased protection over the wrist and forearm
- Made from soft natural rubber latex for improved elasticity and greater comfort
- Reduced protein rating

• IN COMPLIANCE WITH PERSONAL PROTECTIVE EQUIPMENT DIRECTIVE 89/686/EEC, COMPLEX DESIGN CATEGORY; EUROPEAN MEDICAL DEVICE DIRECTIVE 93/42/EEC, CE CLASS 1; TYPE TESTED TO EN 374-2, EN 374-3, EN 388 & EN 420, AND ALSO COMPLIES WITH EN 455 PARTS 1, 2, AND 3.



Synetron® Plus Attributes



	EN 455	Synetron® Plus
Length (mm/in):	270 / 9.4	290 / 11.4
Cuff Thickness (mils)*:	n/a	6.7
Palm Thickness (mils)*:	min. 4.3	9.1
Finger Thickness (mils)*:	min. 4.3	11.8
Force at Break		
Before Aging (N):	min. 9	21
After Aging (N):	min. 6	18
Elasticity/Elongation:		
Before Aging (%):	min. 650	750
After Aging (%):	min. 500	600
Reduced Protein Rating:	max. 200 µg/dm ²	max. 50 µg/g

Order #
 SSYP-INT-S LSYP-INT-L
 MSYP-INT-M XL . . .SYP-INT-XL



EN-374-2,
AQL 1,5



EN-374



Material: Latex
Type: Non-sterile
Color: Natural
Cuff: Extended, beaded
Glove Interior: Powder-free, polymer-coated
Pkg. S-XL: 500 gloves per case, by weight
 10 boxes per case, 50 gloves per box

Usage: For single use only
Shape: Ambidextrous
Size: S - XL
Glove Exterior: Fully textured

⚠ Caution: This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.

*ISO 11193

See for yourself why Microflex® gloves are the first choice of exam glove users worldwide. Contact a Microflex® distributor today to purchase quality Microflex® products. For the name of a distributor in your area, call: **(775) 746-6600**. Fax: (775) 746-6586.

