

# ULTRASENSE™

## POWDER-FREE NITRILE EXAMINATION GLOVES

### • STANDARD EXAMINATION GLOVE


**UltraComfort...** Ultra-soft, next generation nitrile formulation provides the relaxed fit and feel that you would expect from a latex glove.

**UltraDurable...** High quality nitrile offers reliable barrier protection against many hazardous and infectious substances.

**UltraGrip...** Unique manufacturing process delivers an outstanding grip in wet and dry environments.

**UltraValue...** The comfort of latex in an affordable nitrile glove provides the superior value you have been searching for in a synthetic glove.

EN 374-2, EN 374-3 certified



If you are looking for a synthetic glove with the fit and feel of latex, look no further. **UltraSense™** gloves are made from a next generation nitrile formulation to provide comfort, while maintaining the reliable barrier protection of a synthetic glove.

**UltraSense™** gloves are the best overall choice for comfort and value.

EN 388  
(0 0 2 0)



EN 374



EN 374



# ULTRASENSE™

## POWDER-FREE NITRILE EXAMINATION GLOVES

### STANDARD EXAMINATION GLOVE



### GLOVE PROFILE

#### UltraSense™

<b>Length</b>	245 mm / 9.6 in.
<b>Cuff Thickness</b>	2.8 mils
<b>Palm Thickness</b>	3.5 mils
<b>Finger Thickness</b>	4.7 mils
<b>Force at Break</b>	min. 9
<b>Elasticity/Elongation (%)</b> Before Aging (%)	min. 500
<b>Powder Mass</b>	max. 2 mg / glove

### GLOVE ATTRIBUTES

<b>Material:</b>	Nitrile; not made from natural rubber latex
<b>Type:</b>	Non-sterile
<b>Shape:</b>	Ambidextrous
<b>Color:</b>	Custom blue
<b>Sizes:</b>	Extra-small, Small, Medium, Large, Extra-large
<b>Cuff:</b>	Beaded
<b>Glove Interior:</b>	Powder-free
<b>Glove Exterior:</b>	Textured fingers
<b>Usage:</b>	For single use only
<b>Packaging:</b>	100 gloves per dispenser box, by weight; 10 boxes per case

\* ISO11193-1; EN455 I, II, III; EN 374-2, EN 374-3; CE

### Powder-free synthetic examination gloves

See for yourself why **Microflex®** gloves are the first choice of exam glove users nationwide. For samples or more information on any of the fine latex and synthetic products from **Microflex®**, contact your **Microflex®** distributor or for the name of a distributor in your area, call:

**775-746-6600** ext. 1726

Or visit our website at:

[www.microflex.com/int](http://www.microflex.com/int)

Fax: 775-746-6586

PRODUCT CODE	RANGE OF SIZE	DESCRIPTION	CASE QUANTITY
US-INT-XS	EXTRA-SMALL	Nitrile, Non-sterile, Examination	1000
US-INT-S	SMALL	Nitrile, Non-sterile, Examination	1000
US-INT-M	MEDIUM	Nitrile, Non-sterile, Examination	1000
US-INT-L	LARGE	Nitrile, Non-sterile, Examination	1000
US-INT-XL	EXTRA-LARGE	Nitrile, Non-sterile, Examination	1000

**Caution Statement:** Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.



©2007 Microflex Corporation  
P.O. Box 32000  
Reno, NV 89533-2000  
All Rights Reserved.  
07.00003 AL 061907

PATENT PENDING