

Microflex® Supreno Plus™

POWDER-FREE MEDICAL-GRADE NITRILE EXAM GLOVES



A Leader in the EMS Market

Supreno Plus™ nitrile gloves offer durability and additional protection for the wrist and arm, specifically for high-risk applications including the use of chemotherapy drugs. Designed to maintain a comfortable fit while providing reliable hazard protection, Supreno Plus™ gloves are also adhesive-resistant.

Product Attributes: Blue Color • Nitrile • Non-Sterile • Ambidextrous • Textured Fingers
Extended Beaded Cuff • For Single Use Only




Recommended Industries:



DENTAL EMS HEALTHCARE IND. SAFETY LABORATORY

Microflex Product Codes:



AQL 0,65 (LEVEL 3)

S ... SP-INT-S L ... SP-INT-L
M ... SP-INT-M XL ... SP-INT-XL

Packaging:

50 gloves/box
10 boxes/case

GLOVE PROFILE	CUFF LENGTH mm (inches)	CUFF THICKNESS mils (mm)	PALM THICKNESS mils (mm)	FINGER THICKNESS mils (mm)	TENSILE STRENGTH / FORCE AT BREAK		ELASTICITY / ELONGATION		PINHOLE DEFECT RATE (AQL)
					Before Aging	After Aging	Before Aging (%)	After Aging (%)	
ASTM D6319	220 (8.7)	n/a	min. 2.0	min. 2.0	min. 14 MPa	min. 14 MPa	min. 500	min. 500	2.5 Pinhole Defect Rate is FDA minimum Req.
Supreno Plus™	295 (11.6)	4.7 (0.12)	5.9 (0.15)	7.9 (0.2)	30 MPa	27 MPa	500	500	0.65
EN 455	240 (9.45)	n/a	n/a	n/a	6N	6N	n/a	n/a	1.5
Supreno Plus™	295 (11.6)	4.7 (0.12)	5.9 (0.15)	7.9 (0.2)	12N	14N	500	500	0.65

Caution: Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.

The reason is the quality...the difference is Microflex.

See for yourself why Microflex® gloves are the first choice of exam glove users worldwide.
For the name of a distributor in your area, call: **775 746-6600**.
For more information, please visit: www.microflex.com/int

MICROFLEX®
THE MOST TRUSTED NAME IN GLOVES®
A Subsidiary of BarrierSafe Solutions International®, Inc.