

# UltraForm® SE

POWDER-FREE NITRILE EXAMINATION GLOVES

• STANDARD EXAMINATION GLOVE



## Features & Benefits:

- Soft and comfortable for optimal fit and feel
- Forms to your hand to allow full range of motion
- Fully textured for maximum grip in wet or dry conditions
- Non-latex white nitrile
- Meets all current and proposed FDA and EN/ISO standards assuring consistent quality in every box, every time
- EN 374-2, EN 374-3 certified



**Exceptional comfort in a tough nitrile glove.**

**UltraForm®SE** is a soft, comfortable nitrile glove that provides the durability you desire in a non-latex glove. Powder-free **UltraForm®SE** is fully textured for maximum grip in wet or dry conditions.

**UltraForm®SE** works with your hand allowing a full range of motion, while helping to lessen the effect of hand fatigue and stress. **UltraForm®SE** is soft for easy fit and comfort over a long period of time, without compromising integrity and strength.

EN 374



**MICROFLEX®**  
THE MOST TRUSTED NAME IN GLOVES®

[www.microflex.com/int](http://www.microflex.com/int)

REGISTERED  
**ISO 13485**

# UltraForm® SE

POWDER-FREE NITRILE EXAMINATION GLOVES

● STANDARD EXAMINATION GLOVE



**GLOVE PROFILE**

UltraForm®SE\*

|                       |                   |
|-----------------------|-------------------|
| Length                | 240mm / 9.4 in.   |
| Cuff Thickness        | 2.8 mils          |
| Palm Thickness        | 4.3 mils          |
| Finger Thickness      | 5.2 mils          |
| Force at Break        | min. 9            |
| Elasticity/Elongation |                   |
| Before Aging (%):     | min. 500          |
| Powder Mass           | max. 2 mg / glove |

**GLOVE ATTRIBUTES**

|                        |   |
|------------------------|---|
| <b>Material:</b>       | Nitrile; not made from natural rubber latex                   |
| <b>Type:</b>           | Non-sterile   |
| <b>Shape:</b>          | Ambidextrous  |
| <b>Color:</b>          | White   |
| <b>Sizes:</b>          | Extra-Small, Small, Medium, Large, Extra-Large                |
| <b>Cuff:</b>           | Beaded  |
| <b>Glove Interior:</b> | Powder-free   |
| <b>Glove Exterior:</b> | Fully textured  |
| <b>Usage:</b>          | For single use only   |
| <b>Packaging:</b>      | 100 gloves per dispenser box, by weight;<br>10 boxes per case |

\* ISO11193-1; EN455 I, II, III; EN 374-2, EN 374-3; CE Class I

**Powder-free synthetic examination gloves**

See for yourself why **Microflex®** gloves are the first choice of exam glove users nationwide. For samples or more information on any of the fine latex and synthetic products from **Microflex®**, contact your **Microflex®** distributor or for the name of a distributor in your area, call:

**775-746-6600** ext. 1726  
Or visit our website at:  
[www.microflex.com/int](http://www.microflex.com/int)  
Fax: 775-746-6586

| PRODUCT CODE | RANGE OF SIZE | DESCRIPTION                       | CASE QUANTITY |
|--------------|---------------|-----------------------------------|---------------|
| UF-INT-XS    | EXTRA-SMALL   | Nitrile, Non-sterile, Examination | 1000          |
| UF-INT-S     | SMALL         | Nitrile, Non-sterile, Examination | 1000          |
| UF-INT-M     | MEDIUM        | Nitrile, Non-sterile, Examination | 1000          |
| UF-INT-L     | LARGE         | Nitrile, Non-sterile, Examination | 1000          |
| UF-INT-XL    | EXTRA-LARGE   | Nitrile, Non-sterile, Examination | 1000          |

**Caution Statement:** Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.



**MICROFLEX®**  
THE MOST TRUSTED NAME IN GLOVES®

©2007 Microflex Corporation  
P.O.Box 32000  
Reno, NV 89533-2000  
All Rights Reserved.  
07.00003 AL 061907

EUROPEAN PATENT NO. 0456333.