


Microflex® Diamond Grip Plus™

POWDER-FREE MEDICAL-GRADE LATEX EXAM GLOVES



A Leader in Dental, Laboratory and Healthcare Markets


Medical professionals demand the combination of comfort and strength that Diamond Grip Plus™ gloves provide in both wet and dry conditions. Featuring great elasticity and excellent tactile sensitivity, Diamond Grip Plus™ premium multi-purpose gloves deliver continuous and consistent comfort and strength.

Product Attributes: Natural Color • Latex • Non-Sterile • Ambidextrous • Fully Textured Beaded Cuff • For Single Use Only

Recommended Industries:



Microflex Product Codes:



XS ... DGP-INT-XS
S ... DGP-INT-S
M ... DGP-INT-M

L ... DGP-INT-L
XL ... DGP-INT-XL

Packaging:

100 gloves/box
10 boxes/case

GLOVE PROFILE	CUFF LENGTH mm (Inches)	CUFF THICKNESS mils (mm)	PALM THICKNESS mils (mm)	FINGER THICKNESS mils (mm)	TENSILE STRENGTH / FORCE AT BREAK		ELASTICITY / ELONGATION		PINHOLE DEFECT RATE (AQL)
					Before Aging	After Aging	Before Aging (%)	After Aging (%)	
ASTM D3578	220 (8.7)	n/a	min. 3.1	min. 3.1	min. 18 MPa	min. 14 MPa	min. 650	min. 500	2.5 Pinhole Defect Rate is FDA minimum Req.
Diamond Grip Plus™	245 (9.6)	3.9 (0.1)	5.1 (0.13)	6.3 (0.16)	22 MPa	20 MPa	700	600	1.5
EN 455	240 (9.45)	n/a	n/a	n/a	6N	6N	n/a	n/a	1.5
Diamond Grip Plus™	245 (9.6)	3.9 (0.1)	5.1 (0.13)	6.3 (0.16)	9N	6N	700	600	1.5

Caution: This product contains natural rubber latex (latex) which may cause allergic reactions in some users. Safe use of this glove by or on latex sensitized individuals has not been established.

The reason is the quality...the difference is Microflex.

See for yourself why Microflex® gloves are the first choice of exam glove users worldwide. For the name of a distributor in your area, call: **775 746-6600**. For more information, please visit: www.microflex.com/int

MICROFLEX®
THE MOST TRUSTED NAME IN GLOVES®
A Subsidiary of BarrierSafe Solutions International®, Inc.