

Microflex® NeoPro®

POWDER-FREE MEDICAL-GRADE **CHLOROPRENE** EXAM GLOVES



Approved for use with chemotherapy drugs

If your routine calls for secure protection against a variety of chemicals, NeoPro® gloves are your best option. Made with a soft chloroprene blend, NeoPro® gloves offer better puncture resistance and chemical barrier protection than natural rubber latex.

Product Attributes: Green Color • Chloroprene • Non-Sterile • Ambidextrous
Textured Fingers • Beaded Cuff • For Single Use Only



Recommended Industries:



Microflex Product Codes:

XS ... NPG-INT-XS L ... NPG-INT-L
S ... NPG-INT-S XL ... NPG-INT-XL
M ... NPG-INT-M

Packaging:

100 gloves/box
10 boxes/case

GLOVE PROFILE	CUFF LENGTH mm (inches)	CUFF THICKNESS mils (mm)	PALM THICKNESS mils (mm)	FINGER THICKNESS mils (mm)	TENSILE STRENGTH / FORCE AT BREAK		ELASTICITY / ELONGATION		PINHOLE DEFECT RATE (AQL)
					Before Aging	After Aging	Before Aging (%)	After Aging (%)	
ASTM D6319	220 (8.7)	n/a	min. 2.0	min. 2.0	min. 14 MPa	min. 14 MPa	min. 500	min. 400	2.5 <small>Pinhole Defect Rate is FDA minimum Req.</small>
NeoPro®	245 (9.6)	4.3 (0.11)	5.1 (0.13)	6.7 (0.17)	22 MPa	23 MPa	500	400	1.5
EN 455	240 (9.45)	n/a	n/a	n/a	6N	6N	n/a	n/a	1.5
NeoPro®	245 (9.6)	4.3 (0.11)	5.1 (0.13)	6.7 (0.17)	9N	9N	500	400	1.5

Caution: Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.

The reason is the quality...the difference is Microflex.

See for yourself why Microflex® gloves are the first choice of exam glove users worldwide.
For the name of a distributor in your area, call: **775 746-6600**.
For more information, please visit: www.microflex.com/int

MICROFLEX®
THE MOST TRUSTED NAME IN GLOVES®
A Subsidiary of BarrierSafe Solutions International®, Inc.