

Microflex® SafeGrip®

POWDER-FREE EXTENDED CUFF MEDICAL-GRADE LATEX EXAM GLOVES



A Leader in Industrial Safety, Automotive and EMS Markets

For the protection required in high-risk or heavy-duty conditions, rely on SafeGrip® latex gloves. Four times thicker than standard latex gloves, SafeGrip® provides unparalleled defense against many hazards in the work place.

Product Attributes: Blue Color • Latex • Non-Sterile • Ambidextrous • Textured Fingertips
Extended Beaded Cuff • For Single Use Only

Recommended Industries:



Microflex Product Codes:

S ... SG-375-S L SG-375-L
M SG-375-M XL ... SG-375-XL

Packaging:

50 gloves/box
10 boxes/case

GLOVE PROFILE	CUFF LENGTH mm (inches)	CUFF THICKNESS mils (mm)	PALM THICKNESS mils (mm)	FINGER THICKNESS mils (mm)	TENSILE STRENGTH / FORCE AT BREAK		ELASTICITY / ELONGATION		PINHOLE DEFECT RATE (AQL)
					Before Aging	After Aging	Before Aging (%)	After Aging (%)	
ASTM D3578	220 (8.7)	n/a	min. 3.1	min. 3.1	min. 18 MPa	min. 14 MPa	min. 650	min. 500	2.5 <small>Pinhole Defect Rate is FDA minimum Req.</small>
SafeGrip®	295 (11.6)	7.9 (0.20)	11.0 (0.28)	14.2 (0.36)	34 MPa	28 MPa	750	600	1.5
EN 455	240 (9.45)	n/a	n/a	n/a	6N	6N	n/a	n/a	1.5
SafeGrip®	295 (11.6)	7.9 (0.20)	11.0 (0.28)	14.2 (0.36)	30N	15N	750	600	1.5

Caution: This product contains natural rubber latex (latex) which may cause allergic reactions in some users. Safe use of this glove by or on latex sensitized individuals has not been established.

The reason is the quality...the difference is Microflex.

See for yourself why Microflex® gloves are the first choice of exam glove users worldwide. For the name of a distributor in your area, call: **775 746-6600**. For more information, please visit: www.microflex.com/int

MICROFLEX®
THE MOST TRUSTED NAME IN GLOVES®
A Subsidiary of BarrierSafe Solutions International®, Inc.