

# Microflex® Synetron® Plus

POWDER-FREE MEDICAL-GRADE LATEX EXAM GLOVES



## A Leader in the Emergency Medical Services Market


Polymer-coated to slide on easily in wet or dry conditions, Synetron® Plus gloves are fully textured for a secure grip. Exceptional durability to handle tough tasks coupled with fingertip sensitivity for finding a pulse and starting an IV.

**Product Attributes:** Natural Color • Latex • Non-Sterile • Ambidextrous • Fully Textured  
Extended Beaded Cuff • For Single Use Only

### Recommended Industries:



### Microflex Product Codes:



(AQL 1,5) LEVEL 2

S ... SYP-INT-S  
M ... SYP-INT-M

L ... SYP-INT-L  
XL ... SYP-INT-XL

### Packaging:

50 gloves/box  
10 boxes/case

GLOVE PROFILE	CUFF LENGTH mm (inches)	CUFF THICKNESS mils (mm)	PALM THICKNESS mils (mm)	FINGER THICKNESS mils (mm)	TENSILE STRENGTH / FORCE AT BREAK		ELASTICITY / ELONGATION		PINHOLE DEFECT RATE (AQL)
					Before Aging	After Aging	Before Aging (%)	After Aging (%)	
ASTM D3578	220 (8.7)	n/a	min. 3.1	min. 3.1	min. 18 MPa	min. 14 MPa	min. 650	min. 500	2.5 Pinhole Defect Rate is FDA minimum Req.
Synetron® Plus	290 (11.4)	5.9 (0.15)	8.3 (0.21)	10.2 (0.26)	26 MPa	22 MPa	700	600	1.5
EN 455	240 (9.45)	n/a	n/a	n/a	6N	6N	n/a	n/a	1.5
Synetron® Plus	290 (11.4)	5.9 (0.15)	8.3 (0.21)	10.2 (0.26)	15N	11N	700	600	1.5

**Caution:** This product contains natural rubber latex (latex) which may cause allergic reactions in some users. Safe use of this glove by or on latex sensitized individuals has not been established.

**The reason is the quality...the difference is Microflex.**

See for yourself why Microflex® gloves are the first choice of exam glove users worldwide. For the name of a distributor in your area, call: **775 746-6600**. For more information, please visit: [www.microflex.com/int](http://www.microflex.com/int)

**MICROFLEX®**  
THE MOST TRUSTED NAME IN GLOVES®  
A Subsidiary of BarrierSafe Solutions International®, Inc.