

Microflex® UltraSense®

POWDER-FREE MEDICAL-GRADE NITRILE EXAM GLOVES




A Leader in Dental, Laboratory and Healthcare Markets

The non-stick materials in UltraSense® gloves keep you clear of adhesives while still providing a precise grip in both wet and dry conditions. UltraSense® gloves are designed to provide comfort even during extended use with enhanced tactile sensitivity, providing you with a more reliable grip.

Product Attributes: Custom Blue Color • Nitrile • Non-Sterile • Ambidextrous
Textured Fingers • Beaded Cuff • For Single Use Only



Recommended Industries:



(AQL 1,5) LEVEL 2

Microflex Product Codes:

XS ... US-INT-XS L ... US-INT-L
S ... US-INT-S XL ... US-INT-XL
M ... US-INT-M

Packaging:

100 gloves/box
10 boxes/case

GLOVE PROFILE	CUFF LENGTH mm (inches)	CUFF THICKNESS mils (mm)	PALM THICKNESS mils (mm)	FINGER THICKNESS mils (mm)	TENSILE STRENGTH / FORCE AT BREAK		ELASTICITY / ELONGATION		PINHOLE DEFECT RATE (AQL)
					Before Aging	After Aging	Before Aging (%)	After Aging (%)	
ASTM D6319	220 (8.7)	n/a	min. 2.0	min. 2.0	min. 14 MPa	min. 14 MPa	min. 500	min. 400	2.5 Pinhole Defect Rate is FDA minimum Req.
UltraSense®	240 (9.45)	2.8 (0.07)	3.1 (0.08)	4.3 (0.11)	30 MPa	25 MPa	500	400	1.5
EN 455	240 (9.45)	n/a	n/a	n/a	6N	6N	n/a	n/a	1.5
UltraSense®	240 (9.45)	2.8 (0.07)	3.1 (0.08)	4.3 (0.11)	7N	7N	500	400	1.5

Caution: Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.

The reason is the quality...the difference is Microflex.

See for yourself why Microflex® gloves are the first choice of exam glove users worldwide.
For the name of a distributor in your area, call: **775 746-6600**.
For more information, please visit: www.microflex.com/int

MICROFLEX®
THE MOST TRUSTED NAME IN GLOVES®
A Subsidiary of BarrierSafe Solutions International®, Inc.