

# Microflex® UltraSense® Plus

POWDER-FREE MEDICAL-GRADE NITRILE EXAM GLOVES



**Packaged in boxes of 100 for a better value**

UltraSense® Plus gloves provide reliable barrier protection, with enhanced tactile sensitivity, better grip and extended coverage of the wrist and forearm. UltraSense® Plus gloves have a non-stick surface that prevents tapes and adhesives from sticking.


**Product Attributes:** Custom Blue Color • Nitrile • Non-Sterile • Ambidextrous  
Textured Fingers • Extended Beaded Cuff • For Single Use Only



## Recommended Industries:



DENTAL EMS HEALTHCARE IND. SAFETY LABORATORY



(AQL 1,5) LEVEL 2

## Microflex Product Codes:

XS ... USP-INT-XS  
S ... USP-INT-S  
M ... USP-INT-M

L ... USP-INT-L  
XL ... USP-INT-XL

## Packaging:

100 gloves/box  
10 boxes/case

GLOVE PROFILE	CUFF LENGTH mm (Inches)	CUFF THICKNESS mils (mm)	PALM THICKNESS mils (mm)	FINGER THICKNESS mils (mm)	TENSILE STRENGTH / FORCE AT BREAK		ELASTICITY / ELONGATION		PINHOLE DEFECT RATE (AQL)
					Before Aging	After Aging	Before Aging (%)	After Aging (%)	
ASTM D6319	220 (8.7)	n/a	min. 2.0	min. 2.0	min. 14 MPa	min. 14 MPa	min. 500	min. 400	2.5 Pinhole Defect Rate is FDA minimum Req.
UltraSense® Plus	290 (11.4)	2.8 (0.07)	3.5 (0.09)	4.7 (0.12)	33 MPa	32 MPa	500	500	1.5
EN 455	240 (9.45)	n/a	n/a	n/a	6N	6N	n/a	n/a	1.5
UltraSense® Plus	290 (11.4)	2.8 (0.07)	3.5 (0.09)	4.7 (0.12)	7N	7N	500	500	1.5

**Caution:** Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.

**The reason is the quality...the difference is Microflex.**

See for yourself why Microflex® gloves are the first choice of exam glove users worldwide.  
For the name of a distributor in your area, call: **775 746-6600**.  
For more information, please visit: [www.microflex.com/int](http://www.microflex.com/int)

**MICROFLEX®**  
THE MOST TRUSTED NAME IN GLOVES®  
A Subsidiary of BarrierSafe Solutions International®, Inc.