

Microflex®  
**TRANQUILITY®**

**Powder-Free White  
 Nitrile Exam Gloves**

**Item #**

- Extra Small..... TQ-INT-XS
- Small ..... TQ-INT-S
- Medium..... TQ-INT-M
- Large..... TQ-INT-L
- Extra Large..... TQ-INT-XL



**Recommended for use in the following industries:**  
 Dental / Healthcare / Industrial Safety / Laboratory

**Rejuvenate your hands with  
 Tranquility® gloves!**

For the perfect balance of comfort and protection, choose **Tranquility®** white nitrile gloves. **Tranquility®** gloves are latex-free and contain no natural rubber proteins. The glove of choice for individuals that require a soft, thin, comfortable glove for everyday use. **Tranquility®** gloves conform to the hand allowing for excellent tactile sensitivity and grip. The compact box size of **Tranquility®** gloves saves on storage requirements.

**Recommended for use when:**

- A non-latex solution is required, at a value price
- Enhanced comfort and tactile sensitivity are a must
- Everyday use of small instruments demands superior wet and dry grip



**In Compliance with:**

- Personal Protective Equipment 89/686/EEC
- Cat. III, Complex Design
- EN-374-2; Level 2; AQL 1,5
- EN-374-3
- Medical Device Directive 93/42/EEC
- EN-455 Parts 1, 2 and 3
- ASTM F1671

Microflex® Tranquility® Profile		For current specification information, please visit <a href="http://www.microflex.com">www.microflex.com</a> .	
	Tranquility®		
<b>Length (mm):</b>	240	<b>Material:</b> Nitrile	<b>Type:</b> Non-sterile
<b>Cuff Thickness (mm):</b>	0,07	<b>Size:</b> XS - XL	<b>Color:</b> White
<b>Palm Thickness (mm):</b>	0,08	<b>Shape:</b> Ambidextrous	<b>Usage:</b> For single use only
<b>Finger Thickness (mm):</b>	0,11	<b>Cuff:</b> Standard examination length (9.5"), Beaded	
<b>Force at Break:</b>		<b>Glove Exterior:</b> Textured Fingers	
Before Aging (Newtons):	7	<b>Glove Interior:</b> Powder-Free, Chlorinated	
After Aging (Newtons):	7	<b>Packaging:</b> 1,000 gloves/case, by weight; 10 boxes/case	
<b>Elasticity/Elongation:</b>		100 gloves per dispenser box, by weight	
Before Aging (%):	500	<b>Caution:</b> Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.	
After Aging (%):	400		

**The reason is the quality... the difference is Microflex®.**  
 See for yourself why Microflex® gloves are the first choice of exam glove users nationwide. For **FREE SAMPLES**, please visit [www.international.microflex.com](http://www.international.microflex.com). For the name of a distributor in your area, call: **(775) 746-6600**.

