

Microflex®

TRANQUILITY®

Powder-Free White Nitrile* Exam Gloves

Item

Extra Small.....	TQ-601-XS
Small.....	TQ-601-S
Medium.....	TQ-601-M
Large.....	TQ-601-L
Extra Large.....	TQ-601-XL



Recommend for use in the following industries:

Healthcare / Dental / Industrial Safety / Laboratory

Rejuvenate your hands with Tranquility® gloves!

For the perfect balance of comfort and protection, choose **Tranquility®** white nitrile. **Tranquility®** gloves are latex-free and contain no natural rubber proteins. The glove of choice for individuals that require a soft, thin, comfortable glove for everyday use. **Tranquility®** gloves conform to the hand allowing for excellent tactile sensitivity and grip. **Tranquility's** compact box size saves on storage requirements.

Recommended for use when:

- A latex-free solution is required, at a value price
- Enhanced comfort and tactile sensitivity are a must
- Everyday use of small instruments demands superior wet and dry grip

NOT MADE FROM NATURAL RUBBER LATEX

*U.S. Patent No. 6,451,893, U.S. Patent No. 7,176,260

Microflex® Tranquility® Profile			For current specification information, please visit www.microflex.com .	
	ASTM D 6319	Tranquility®		
Length (mm/in):	220 / 8.7	240 / 9.6	Material: Nitrile	Type: Non-sterile
Cuff Thickness (mils):	n/a	2.8 ± 0.8	Size: XS - XL	Color: White
Palm Thickness (mils):	min. 2.0	3.1 ± 0.5	Shape: Ambidextrous	Usage: For single use only
Finger Thickness (mils):	min. 2.0	4.3 ± 1.2	Cuff: Standard examination length (9.5"), Beaded	
Tensile Strength:			Glove Exterior: Textured fingers	
Before Aging (MPa):	min. 14	min. 18	Glove Interior: Powder-free, chlorinated	
After Aging (MPa):	min. 14	min. 16	Packaging: 1,000 gloves/case, by weight; 10 boxes/case	
Elasticity/Elongation:			100 gloves per dispenser box, by weight	
Before Aging (%):	min. 500	min. 500	Caution: Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.	
After Aging (%):	min. 400	min. 400		

The reason is the quality... the difference is Microflex®.

See for yourself why Microflex® gloves are the first choice of exam glove users nationwide. For **FREE SAMPLES**, please visit www.microflex.com. For the name of distributor in your area, please **CALL (1-800-876-6866)** 6am-5pm (pst) or **FAX (1-800-876-6632)** 24 hours

MICROFLEX®
THE MOST TRUSTED NAME IN GLOVES®
A Division of BarrierSafe Solutions International, Inc.