

Microflex®
DIAMOND GRIP PLUS™

Powder-Free **Latex**
 Exam Gloves

Item #

- Extra Small..... DGP-INT-XS
- Small DGP-INT-S
- Medium..... DGP-INT-M
- Large..... DGP-INT-L
- Extra Large..... DGP-INT-XL



Recommended for use in the following industries:

Acute & Alternate Care / Automotive / Dental / EMS / Foodservice & Processing / Industrial Safety / Laboratory / Law & Corrections

The combination of comfort and strength make **Diamond Grip Plus™** a leading latex glove of choice for many medical professionals. **Diamond Grip Plus™** gloves offer a consistent and secure wet and dry grip. Great elasticity, a comfortable fit, and excellent tactile sensitivity make **Diamond Grip Plus™** a perfect multi-purpose glove.

Diamond Grip Plus™ gloves offer a superior combination of comfort and durability in an all-purpose glove.

Recommended for use when:

- Durable comfort is required for lasting wear
- Consistent grip is necessary for handling utensils on a day-to-day basis
- Tactile sensitivity is important when working with a number of different instruments

EN374-2: 2003 EN374-3: 2003



LEVEL 2

In Compliance with:

- Personal Protective Equipment 89/686/EEC
- Cat. III, Complex Design
- EN-374-2; Level 2; AQL 1,5
- EN-374-3
- Medical Device Directive 93/42/EEC
- EN-455 Parts 1, 2 and 3
- ASTM F1671

Microflex® Diamond Grip Plus™ Profile		For current specification information, please visit www.microflex.com .	
	Diamond Grip Plus™		
Length (mm):	245	Material: Latex	Type: Non-sterile
Cuff Thickness (mm):	0,10	Size: XS - XL	Color: Natural
Palm Thickness (mm):	0,13	Shape: Ambidextrous	Usage: For single use only
Finger Thickness (mm):	0,16	Cuff: Standard examination length (9.5"), Beaded	
Force at Break:		Glove Exterior: Fully Textured	
Before Aging (Newtons):	9	Glove Interior: Powder-Free, Chlorinated	
After Aging (Newtons):	6	Packaging: 1,000 gloves/case, by weight; 10 boxes/case	
Elasticity/Elongation:		100 gloves per dispenser box, by weight	
Before Aging (%):	700		
After Aging (%):	600		
		Caution: This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.	

The reason is the quality... the difference is Microflex®. See for yourself why Microflex® gloves are the first choice of exam glove users nationwide. For **FREE SAMPLES**, please visit www.international.microflex.com. For the name of a distributor in your area, call: **(775) 746-6600**.

