

Microflex® FREEFORM® SE

Powder-Free Blue Nitrile Exam Gloves

Item

Extra Small.....	FFS-INT-XS
Small	FFS-INT-S
Medium.....	FFS-INT-M
Large.....	FFS-INT-L
Extra Large.....	FFS-INT-XL



Recommended for use in the following industries:

Acute & Alternate Care / Automotive / Dental / EMS / Foodservice & Processing / Industrial Safety / Laboratory / Law & Corrections

If your need is softness combined with durability, **FreeForm® SE** is the glove for you. Designed with strength and comfort in mind, **FreeForm® SE** gloves are a unique combination of features in a nitrile glove.

FreeForm® SE gloves are the best overall value in a comfortable mid-range synthetic glove.

Recommended for use when:

- Working in wet environments that demand a reliable grip
- Light to mid-duty tasks require a combination of comfort and durability
- Precise contact and softness require enhanced tactile sensitivity

Features and benefits:

- Non-sitck surface properties resistant to tapes and adhesives

EN374-2: 2003 EN374-3: 2003



LEVEL 2



In Compliance with:

- Personal Protective Equipment 89/686/EEC
- Cat. III, Complex Design
- EN-374-2; Level 2; AQL 1,5
- EN-374-3
- Medical Device Directive 93/42/EEC
- EN-455 Parts 1, 2 and 3
- ASTM F1671

Microflex® FreeForm® SE Profile		For current specification information, please visit www.microflex.com .	
	FreeForm® SE		
Length (mm):	245	Material: Nitrile	Type: Non-sterile
Cuff Thickness (mm):	0,07	Size: XS - XL	Color: Blue
Palm Thickness (mm):	0,09	Shape: Ambidextrous	Usage: For single use only
Finger Thickness (mm):	0,13	Cuff: Standard examination length (9.5"), Beaded	
Force at Break:		Glove Exterior: Textured Fingers, Polymer-Coated	
Before Aging (Newtons):	10	Glove Interior: Powder-Free, Chlorinated	
After Aging (Newtons):	10	Packaging: 1,000 gloves/case, by weight; 10 boxes/case	
Elasticity/Elongation:		100 gloves per dispenser box, by weight	
Before Aging (%):	500	Caution: Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.	
After Aging (%):	400		

The reason is the quality... the difference is Microflex®. See for yourself why Microflex® gloves are the first choice of exam glove users nationwide. For **FREE SAMPLES**, please visit www.international.microflex.com. For the name of a distributor in your area, call: **(775) 746-6600**.

MICROFLEX®
THE MOST TRUSTED NAME IN GLOVES®
A division of BarrierSafe Solutions International, Inc.