

Microflex®

# MICRO ONE®

## Lightly Powdered Latex Exam Gloves

### Item #

Extra Small.....	MO-INT-XS
Small .....	MO-INT-S
Medium.....	MO-INT-M
Large.....	MO-INT-L
Extra Large.....	MO-INT-XL



### Recommended for use in the following industries:

Acute & Alternate Care / Automotive / Dental / EMS / Foodservice & Processing / Industrial Safety / Laboratory / Law & Corrections

The **Micro One**® glove is a highly versatile examination glove that performs well in demanding conditions. Made from soft latex, **Micro One**® gloves have a natural fit for reduced hand fatigue. Lightly powdered **Micro One**® gloves don easily and have a smooth surface for uncompromised sensitivity. **Micro One**® gloves are thicker and longer than typical latex exam gloves for added protection where you need it most. **Micro One**® gloves meet or exceed current EN standards for exam gloves quality and performance.

### Features and benefits:

- Longer than typical latex examination gloves for added protection where you need it most
- Smooth surface for excellent tactile sensitivity
- Meets maximum powder and maximum water extractable protein content guidelines sensitivity

### In Compliance with:

- Medical Device Directive 93/42/EEC
- EN-455 Parts 1, 2 and 3
- ASTM F1671

### Microflex® Micro One® Profile

For current specification information, please visit [www.microflex.com](http://www.microflex.com).

Micro One®			
<b>Length (mm):</b>	245	<b>Material:</b> Latex	<b>Type:</b> Non-sterile
<b>Cuff Thickness (mm):</b>	0,10	<b>Size:</b> XS - XL	<b>Color:</b> White
<b>Palm Thickness (mm):</b>	0,14	<b>Shape:</b> Ambidextrous	<b>Usage:</b> For single use only
<b>Finger Thickness (mm):</b>	0,16	<b>Cuff:</b> Standard examination length (9.5"), Beaded	
<b>Force at Break:</b>		<b>Glove Exterior:</b> Smooth	
Before Aging (Newtons):	9	<b>Glove Interior:</b> Lightly Powdered, absorbable cornstarch (USP)	
After Aging (Newtons):	6	<b>Packaging:</b> 1,000 gloves/case, by weight; 10 boxes/case	
<b>Elasticity/Elongation:</b>		100 gloves per dispenser box, by weight	
Before Aging (%):	750	<b>Caution:</b> This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.	
After Aging (%):	600		

### The reason is the quality... the difference is Microflex®.

See for yourself why Microflex® gloves are the first choice of exam glove users nationwide. For **FREE SAMPLES**, please visit [www.international.microflex.com](http://www.international.microflex.com). For the name of a distributor in your area, call: **(775) 746-6600**.

**MICROFLEX®**  
THE MOST TRUSTED NAME IN GLOVES®  
A division of BarrierSafe Solutions International, Inc.