

Microflex®

# SUPRENO PLUS™

## Powder-Free Blue-Violet Nitrile Exam Gloves

### Item #

- Extra Small..... SP-INT-XS
- Small ..... SP-INT-S
- Medium..... SP-INT-M
- Large..... SP-INT-L
- Extra Large..... SP-INT-XL



**Recommended for use in the following industries:**  
Automotive / Dental / EMS / Healthcare / Industrial Safety / Laboratory

If your need is strength and durability with extra protection, **Supreno Plus™** gloves are an excellent choice. With enhanced puncture resistance and exceptional toughness, latex-free **Supreno Plus™** gloves provide reliable protection for high-risk applications while maintaining a comfortable fit. **Supreno Plus™** gloves offer the protection you need against the toughest hazards you face on the job, day in and day out.

With an extended cuff for added protection over the wrist and arm, **Supreno Plus™** gloves are the best combination of strength and durability in an allpurpose non-latex glove.

**Recommended for use when:**

- Contact with specific chemicals is a standard activity each day\*
- Strength and protection are of great importance
- Day-to-day hazards require the best combination of durability and security

\* Tested for use with chemotherapy drugs. Gloves used for protection against chemotherapy drugs exposure must be selected specifically for the type of chemicals used. Tested for use against Carmustine, Cyclophosphamide, Doxorubicin HCl, 5-Fluorouracil, Cisplatin, Etoposide, Paclitaxel, Thio-Tepa, Mitomycin C, Daunorubicin HCl, Ifosfamide, Melphalan, Idarubicin HCl, Carboplatin, Bleomycin Sulfate, Mitoxantrone, Cytarabine HCl



### In Compliance with:

- Personal Protective Equipment 89/686/EEC
- Cat. III, Complex Design
- EN-374-2; Level 3; AQL 0,65
- EN-374-3
- Medical Device Directive 93/42/EEC
- EN-455 Parts 1, 2 and 3
- ASTM F1671

Microflex® Supreno Plus™ Profile		For current specification information, please visit <a href="http://www.microflex.com">www.microflex.com</a> .	
	Supreno Plus™		
<b>Length (mm):</b>	295	<b>Material:</b> Nitrile	<b>Type:</b> Non-sterile
<b>Cuff Thickness (mm):</b>	0,12	<b>Size:</b> S - XXL	<b>Color:</b> Blue-Violet
<b>Palm Thickness (mm):</b>	0,15	<b>Shape:</b> Ambidextrous	<b>Usage:</b> For single use only
<b>Finger Thickness (mm):</b>	0,20	<b>Cuff:</b> Standard examination length (9.5"), Beaded	
<b>Force at Break:</b>		<b>Glove Exterior:</b> Textured Fingers, Polymer-Coated	
Before Aging (Newtons):	12	<b>Glove Interior:</b> Powder-Free, Chlorinated	
After Aging (Newtons):	14	<b>Packaging:</b> 1,000 gloves/case, by weight; 10 boxes/case	
<b>Elasticity/Elongation:</b>		100 gloves per dispenser box, by weight	
Before Aging (%):	500		
After Aging (%):	400		
<b>Caution:</b> Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.			

**The reason is the quality... the difference is Microflex®.** See for yourself why Microflex® gloves are the first choice of exam glove users nationwide. For **FREE SAMPLES**, please visit [www.international.microflex.com](http://www.international.microflex.com). For the name of a distributor in your area, call: **(775) 746-6600**.

