

Microflex®

SYNETRON® PLUS

Powder-Free Latex Exam Gloves

Item

- Small SYP-INT-S
- Medium SYP-INT-M
- Large SYP-INT-L
- Extra Large SYP-INT-XL



Recommended for use in the following industries:

EMS / Industrial Safety / Laboratory

Synetron® Plus latex examination gloves are a strong barrier against many chemicals and biohazards, and its comfortable during extended periods of wear. **Synetron® Plus** is textured for a secure grip and is polymer-coated for product consistency and easy donning.

Features and benefits:

- Interior/exterior polymer-coating system provides superior anti-sticking properties that assist with highly effective donning in wet, damp or dry conditions
- Extended cuff for increased protection over the wrist and forearm
- Made from soft natural rubber latex for improved elasticity and greater comfortprotein content guidelines sensitivity
- Reduced protein rating

In Compliance with:

- Medical Device Directive 93/42/EEC
- EN-455 Parts 1, 2 and 3
- ASTM F1671

Microflex® Synetron® Plus Profile

For current specification information, please visit www.microflex.com.

Synetron® Plus			
Length (mm):	290	Material: Latex	Type: Non-sterile
Cuff Thickness (mm):	0,11	Size: S - XL	Color: Natural
Palm Thickness (mm):	0,14	Shape: Ambidextrous	Usage: For single use only
Finger Thickness (mm):	0,18	Cuff: Extended Cuff length (12"), Beaded	
Force at Break:		Glove Exterior: Fully Textured, Polymer-Coated	
Before Aging (Newtons):	15	Glove Interior: Powder-Free, Chlorinated	
After Aging (Newtons):	11	Packaging: 500 gloves/case, by weight; 10 boxes/case	
Elasticity/Elongation:		50 gloves per dispenser box, by weight	
Before Aging (%):	700	Caution: This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.	
After Aging (%):	600		

The reason is the quality... the difference is Microflex®. See for yourself why Microflex® gloves are the first choice of exam glove users nationwide. For **FREE SAMPLES**, please visit www.international.microflex.com. For the name of a distributor in your area, call: **(775) 746-6600**.

MICROFLEX®
THE MOST TRUSTED NAME IN GLOVES®
A division of BarrierSafe Solutions International, Inc.