

Microflex®

UltraSense®

Powder-Free Blue Nitrile Exam Gloves

Item

- Extra Small..... US-220-XS
- Small US-220-S
- Medium..... US-220-M
- Large..... US-220-L
- Extra Large..... US-220-XL



Recommended for use in the following industries:

Dental / EMS / Foodservice / Healthcare /
Industrial Safety / Laboratory

A nitrile glove with the fit and feel of latex.

UltraSense® gloves are made from a patented next-generation nitrile formulation to provide comfort, while maintaining the reliable barrier protection of a synthetic glove. UltraSense® gloves provides a precise grip in both wet and dry conditions, while offering comfort for extended periods of use.

Recommended for use when:

- Handling small instruments, where fit and feel is key
- Light tasks require excellent grip on a day-to-day basis
- Performing extended or repetitive duties that may cause hand fatigue
- When your day-to-day duties require a glove that is resistant to tapes and adhesives

In Compliance with:

- Personal Protective Equipment 89/686/EEC
- Cat. III, Complex Design
- EN-374-2; Level 2; AQL 1,5
- EN-374-3
- Medical Device Directive 93/42/EEC
- EN-455 Parts 1, 2 and 3
- ASTM F1671

Microflex® UltraSense® Profile		For current specification information, please visit www.microflex.com .	
	UltraSense®		
Length (mm):	243	Material: Nitrile	Type: Non-sterile
Cuff Thickness (mm):	0,07	Size: XS - XL	Color: Blue
Palm Thickness (mm):	0,08	Shape: Ambidextrous	Usage: For single use only
Finger Thickness (mm):	0,11	Cuff: Standard examination length (9.5"), Beaded	
Force at Break:		Glove Exterior: Textured Fingers, Polymer-Coated	
Before Aging (Newtons):	7	Glove Interior: Powder-Free, Chlorinated	
After Aging (Newtons):	7	Packaging: 1,000 gloves/case, by weight; 10 boxes/case	
Elasticity/Elongation:		100 gloves per dispenser box, by weight	
Before Aging (%):	500	Caution: Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.	
After Aging (%):	400		

The reason is the quality... the difference is Microflex®. See for yourself why Microflex® gloves are the first choice of exam glove users nationwide. For **FREE SAMPLES**, please visit www.international.microflex.com. For the name of a distributor in your area, call: **(775) 746-6600**.

MICROFLEX®
THE MOST TRUSTED NAME IN GLOVES®
A division of BarrierSafe Solutions International, Inc.