

Microflex®

UltraSense® PLUS

Powder-Free Blue Nitrile
Extended Cuff Exam Gloves

Item

Extra Small..... USP-INT-XS
Small USP-INT-S
Medium..... USP-INT-M
Large..... USP-INT-L
Extra Large..... USP-INT-XL



Recommended for use in the following industries:

Dental / EMS / Laboratory

UltraSense® Plus is a patented, next generation nitrile which provides a reliable barrier while offering the comfort you are used to getting with a latex glove. Engineered with enhanced tactile sensitivity – providing an excellent grip in both wet and dry environments. **UltraSense® Plus** is an extended cuff offering that is packed 100 gloves per box as opposed to the standard 50, maximizing your savings.

Recommended for use when:

- Strength and protection are of great importance
- Working in environments that demand consistent grip and reliability
- Precise contact requires enhanced tactile sensitivity
- When your day-to-day duties require a glove that is resistant to tapes and adhesives

In Compliance with:

- Personal Protective Equipment 89/686/EEC
- Cat. III, Complex Design
- EN-374-2; Level 2; AQL 1,5
- EN-374-3
- Medical Device Directive 93/42/EEC
- EN-455 Parts 1, 2 and 3
- ASTM F1671

Microflex® UltraSense® Plus Profile

For current specification information, please visit www.microflex.com.

UltraSense® Plus			
Length (mm):	290	Material: Nitrile	Type: Non-sterile
Cuff Thickness (mm):	0,07	Size: S - XXL	Color: Blue
Palm Thickness (mm):	0,08	Shape: Ambidextrous	Usage: For single use only
Finger Thickness (mm):	0,11	Cuff: Extended Cuff length (12"), Beaded	
Force at Break:		Glove Exterior: Textured Fingers, Polymer-Coated	
Before Aging (Newtons):	7	Glove Interior: Powder-Free, Chlorinated	
After Aging (Newtons):	7	Packaging: 1,000 gloves/case, by weight; 10 boxes/case	
Elasticity/Elongation:		100 gloves per dispenser box, by weight	
Before Aging (%):	500	Caution: Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.	
After Aging (%):	400		

The reason is the quality... the difference is Microflex®. See for yourself why Microflex® gloves are the first choice of exam glove users nationwide. For **FREE SAMPLES**, please visit www.international.microflex.com. For the name of a distributor in your area, call: **(775) 746-6600**.

MICROFLEX®
THE MOST TRUSTED NAME IN GLOVES®
A division of BarrierSafe Solutions International, Inc.