

Microflex®

# SUPRENO®

## Powder-Free Blue-Violet Nitrile Exam Gloves

### Item #

Extra Small.....	SU-INT-XS
Small .....	SU-INT-S
Medium.....	SU-INT-M
Large.....	SU-INT-L
Extra Large.....	SU-INT-XL



### Recommended for use in the following industries:

Automotive / Dental / EMS / Healthcare / Industrial Safety / Laboratory

If your need is strength and durability, **Supreno®** gloves provide exceptional toughness in a comfortable nitrile glove. **Supreno®** gloves provide reliable protection against the hazards you work with day in and day out.

**Supreno®** gloves are the best combination of strength and durability in an all-purpose synthetic glove.

### Recommended for use when:

- Strength and protection are of great importance
- Day-to-day hazards require the best combination of durability and security
- Working in environments that demand consistent grip and reliability

### In Compliance with:

- Personal Protective Equipment 89/686/EEC
- Cat. III, Complex Design
- EN-374-2; Level 3; AQL 0,65
- EN-374-3
- Medical Device Directive 93/42/EEC
- EN-455 Parts 1, 2 and 3
- ASTM F1671

### Microflex® Supreno® Profile

For current specification information, please visit [www.microflex.com](http://www.microflex.com).

Supreno®			
<b>Length (mm):</b>	245	<b>Material:</b> Nitrile	<b>Type:</b> Non-sterile
<b>Cuff Thickness (mm):</b>	0,10	<b>Size:</b> XS - XL	<b>Color:</b> Blue-Violet
<b>Palm Thickness (mm):</b>	0,12	<b>Shape:</b> Ambidextrous	<b>Usage:</b> For single use only
<b>Finger Thickness (mm):</b>	0,18	<b>Cuff:</b> Standard examination length (9.5"), Beaded	
<b>Force at Break:</b>		<b>Glove Exterior:</b> Textured Fingers	
Before Aging (Newtons):	10	<b>Glove Interior:</b> Powder-Free, Chlorinated	
After Aging (Newtons):	11	<b>Packaging:</b> 1,000 gloves/case, by weight; 10 boxes/case	
<b>Elasticity/Elongation:</b>		100 gloves per dispenser box, by weight	
Before Aging (%):	500	<b>Caution:</b> Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.	
After Aging (%):	400		

### The reason is the quality... the difference is Microflex®.

See for yourself why Microflex® gloves are the first choice of exam glove users nationwide. For **FREE SAMPLES**, please visit [www.international.microflex.com](http://www.international.microflex.com). For the name of a distributor in your area, call: **(775) 746-6600**.

**MICROFLEX®**  
THE MOST TRUSTED NAME IN GLOVES®  
A division of BarrierSafe Solutions International, Inc.